

Researchers worldwide protect more than two billion samples inside Thermo Scientific cold storage equipment. With proven solutions including +4°C refrigerators to -196°C cryogenic freezers you're free to concentrate on your work without worrying about your valuable samples.

Thermo Scientific Harris HD Ultra-Low Temperature Freezers

Store up to 70,000 samples



Easy-to-read microprocessor control delivers industry-leading performance to protect the integrity of your work.



Our HD technology and design allow 31.4 cu. ft. of capacity in a 25 cu. ft. footprint. Store up to 70,000 samples in 35 inventory storage racks.

Features and Benefits

- Up to 20% greater heat removal capacity over ordinary freezers, ensuring faster temperature recovery after door openings*
- Consumes up to 15% less power than other freezers, reducing ownership costs and benefitting the environment*
- Emits up to 15% less heat into the environment than other freezers*
- Noise levels are reduced by up to 40% over previous generation, allowing unit to reside directly in the lab
- Highly-efficient, innovative thin wall insulation technology allows for expanded interior capacity without increasing the exterior footprint
- Microprocessor control with set point controls for temperature and alarm system
- PLUS time delay feature – delays restart of a bank of freezers after building power failure, reducing stress on compressors and electrical components
- On-Board diagnostics report error codes allow the researcher to diagnose problems faster and minimize sample exposure
- Automatic voltage buck/boost system adds to the longevity of the freezer by stabilizing the voltage to the unit
- Optional integrated back-up system provides assurance that samples will not be compromised in event of power outage
- Recessed power switch controls down time due to accidental power shut down
- Reinforced shelves for high-load storage
- Door-mounted pressure equalization port (PEP) does not take up space and does not require maintenance (ice removal)
- Voltage indicator provides fast, accurate information on incoming line voltage
- Optional new communication outputs provide flexibility for connecting to alarm systems
- Battery readout provides user with assurance that alarms and displays will function in the event of a power outage
- Seven-day, 6" circular chart for validation and regulatory requirements

* Based on internal performance data

Thermo Scientific Harris HD Freezers Specifications – Not Available Outside North America

Temperature Range -50°C to -86°C						
Model No.	Capacity cu. ft. (L)	Electrical	Amps/Breaker (Plug)	Interior Dimensions H x W x D inches (mm)	Exterior Dimensions H x W x D inches (mm)	Shipping Weight kg (lbs)
DLT-25V-85HD-A	24.4 (691)	115V/60Hz 1Ø	16/20 (P2)	50.5 x 34.5 x 24.3	77.8 x 43.0 x 36.8	825 (358)
DLT-25V-85HD-D		208-230V/60Hz 1Ø	12/15 (P4)	(1283 x 879 x 616)	(1975 x 1092 x 935)	
DLT-32V-80HD-D	31.4 (889)	208-230V/60Hz 1Ø	12/15 (P4)	50.5 x 39.1 x 27.3	77.8 x 48.0 x 36.8	930 (478)
				(1283 x 993 x 693)	(1975 x 1219 x 935)	

Thermo Scientific Harris Rack Systems

Vials						
Model No.	Description H x F-B x W (inches)	Freezer Vol (cu. ft.)	Qty of 2" Boxes	Qty of 3" Boxes	Racks/Shelf	Racks/Freezer
6113-1	Adjustable Side-Access Rack (9.39 x 22.19 x 5.34)	24.4	16	12	6	30
5992	Sliding Drawer Rack (5.4 x 21.25 x 7)	24.4	12	–	6	30
5996	Sliding Drawer Rack (5.4 x 21.25 x 9.4)	24.4	16	–	6	30
5991	Sliding Drawer Rack (5.4 x 21.25 x 7)	24.4	–	8	6	30
7145	Fixed Side-Access Rack for 2" Boxes (5.4 x 21.25 x 7)	24.4	16	–	6	30
7187	Adjustable Side-Access Rack (9.4 x 26.75 x 5.3)	31.4	20	15	7	35
7243	Sliding Drawer Rack (9.63 x 27.25 x 5.5)	31.4	20	–	7	35

Microplates							
Model No.	Description H x F-B x W (inches)	Freezer Vol (cu. ft.)	Deepwell Plates	Shallow-Well Plates	Microtube Boxes	Racks/Shelf	Racks/Freezer
6733	Microplates, Individual Compartments (8.75 x 22.19 x 5.5)	24.4	66	–	–	6	30
7146	Microplates, Stacked (9.18 x 22 x 5.5)	24.4	84	102	–	6	30
7147	Microplates, Locking Rod (9.37 x 22 x 5.5)	24.4	84	96	–	6	30
7149	Deepwell, Individual Compartments With Locking Rod (9.37 x 22 x 5.5)	24.4	30	–	–	6	30
7150	Microtube Box and Microplate Rack With Lids (8.87 x 22 x 5.5)	24.4	72	–	24	6	30
7152	Microtube Box and Microplate Rack With Lids, Locking Rod (9.06 x 22 x 5.5)	24.4	72	–	24	6	30
7265	Fixed Side-Access Rack for 15-50mL Tubes (5.75 x 5.75 x 4.75)	24.4	–	–	8	6	30
7282	Microplates, Stacked (9.25 x 22.63 x 5.5)	31.4	98	119	–	7	35

Thermo Scientific Harris Inventory Accessories

Accessories	
Model No.	Description
5954	2" High Fiberboard Boxes – 12 qty.
5956	3" High Fiberboard Boxes – 12 qty.
5958	10 x 10, 7/16" Cell – Holds 100 12mm Vials
5959	7 x 7, 5/8" Cell – Holds 49 16mm Vials
5960	8 x 8, 9/16" Cell – Holds 64 14mm Vials
6212	9 x 9, 1/4" Cell – Holds 81 13mm Vials

© 2009 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.